



SOLDER ANALYSIS

315 Fairbank St., Addison, IL 60101

Remove 150 grams or 5 oz. of solder for sampling, directly from the wave in a soldering machine or from a pot after mixing well.

Insert this form into the attached clear envelope.

MAKE SURE THE QUALITEK ADDRESS (OTHER SIDE) IS SHOWING.

Place solder sample in the envelope, staple it closed, affix postage and mail.

Date: _____ P.O.#: _____

Company: _____

Street: _____

City, State, Zip: _____

Results To: _____

Tel. No.: _____ Your Rep./Distributor: _____

E-mail: _____

Visit our website for all alloy options at qualitek.com or call 630-628-8083

MAILER NO.



SOLDER ANALYSIS

315 Fairbank St., Addison, IL 60101

Remove 150 grams or 5 oz. of solder for sampling, directly from the wave in a soldering machine or from a pot after mixing well.

Insert this form into the attached clear envelope.

MAKE SURE THE QUALITEK ADDRESS (OTHER SIDE) IS SHOWING.

Place solder sample in the envelope, staple it closed, affix postage and mail.

Date: _____ P.O.#: _____

Company: _____

Street: _____

City, State, Zip: _____

Results To: _____

Tel. No.: _____ Your Rep./Distributor: _____

E-mail: _____

Visit our website for all alloy options at qualitek.com or call 630-628-8083

MAILER NO.

CUT OFF GREY AREA

FILL OUT CARD, FOLD ALONG DOTTED LINE, PUT IN ENVELOPE AND MAIL

To: **QUALITEK INTERNATIONAL, INC.**

Attn: Quality Assurance Test Lab
315 Fairbank Street
Addison, IL 60101
U.S.A.

CUT HERE FOR TWO ANALYSIS

To: **QUALITEK INTERNATIONAL, INC.**

Attn: Quality Assurance Test Lab
315 Fairbank Street
Addison, IL 60101
U.S.A.

Circle LEADED Option: A B C Alloy: _____

Your Machine # _____ Sample Identifier: _____

Circle LEAD FREE Option: LF-1 LF-2 LF-3 Alloy: _____

Your Machine # _____ Sample Identifier: _____

Qualitek International, Inc. Tel: 630-628-8083 - www.qualitek.com

Circle LEADED Option: A B C Alloy: _____

Your Machine # _____ Sample Identifier: _____

Circle LEAD FREE Option: LF-1 LF-2 LF-3 Alloy: _____

Your Machine # _____ Sample Identifier: _____

Qualitek International, Inc. Tel: 630-628-8083 - www.qualitek.com